



ALMAG-MEDICAL DEVICE

Mahri Satlykova Sahymuradovna

student of the Oguzhan Engineering Technologies University of Turkmenistan
Ashgabat, Turkmenistan

When they reach middle age, most people begin to experience back pain. This problem is associated with a sedentary lifestyle, which today, unfortunately, leads the vast majority of our compatriots. A sitting posture is unusual for a healthy active person and significantly affects the condition of his back, as well as internal organs. Staying in this position for a long time puts the back muscles under excessive strain.

Other causes of back pain can be muscle loads during sports activities, hard work in the country or lifting weights. The consequences are deplorable: unbearable back pain appears, which signals that there have been violations in the musculoskeletal system.

This leads to various diseases: arthritis, osteoarthritis of various joints. These ailments are the worst amenable to drug treatment, and ointments and pills just muffle the pain for a while.

Recently, a new unique device has appeared on sale that can save a person from many diseases. This is "Almag-01" - a device whose action is based on magnetic radiation. The device is very easy to use. The device produces traveling magnetic fluxes that have a beneficial effect on damaged organs.

The use of the device allows you to get rid of back pain, relieve symptoms of inflammation, increase efficiency. Due to its unique design, the device covers a large area of impact: almost the entire spine, upper or lower limbs, torso.

It is quite easy to use Almag-01: the instruction says that the device can be used independently at home. Let's get acquainted with this amazing device.

The device is a small chain of interconnected coils (nodes). Its weight is only 620 grams. One pulse is produced with a duration of 1.5 - 2.5 ms and a frequency of 6 Hz. The device is powered by the mains, when it is turned on, the light comes on. signal.

During application, it is necessary to adhere to the operating mode, which is calculated for the Almag-01 device. The scheme is simple, but it must be followed. The device can work for 6 hours, but every 20 minutes you need to pause. This is necessary in order to avoid overvoltage. The break is maintained for 10 minutes.

The device is intended for treatment:

problems in the musculoskeletal system: osteochondrosis of various localization, arthrosis, bursitis;

damage to the joints or bone tissue of a different nature: injuries, bruises, bruises, sprains or muscles, swelling, burns;

diseases of the cardiovascular system: atherosclerosis of the vessels of the lower extremities, heart murmurs, hypertension;

digestive system problems: gastritis, pancreatitis, ulcers, colitis, and chronic flatulence;

gynecological diseases in women: inflammation, complications after caesarean section, painful menstruation;

venous diseases: thrombosis of the vessels of the lower extremities (chronic and acute), thrombophlebitis, varicose veins;

complications of diabetes mellitus and concomitant diseases;

dermatological problems: healing of scars in the postoperative period, itchy dermatoses;

respiratory diseases: bronchitis, pneumonia, bronchial asthma and other lung ailments;

neurological diseases: vascular dysfunctions of the brain, various disorders of the nervous system.

Magnetic pulses penetrate the tissues of the body to a depth of 8 cm, which allows the device to subtly affect the internal organs, quickly relieve pain and successfully treat various ailments.

The device is rather unusual in shape, so a lot of questions may arise during the first use. What is the best position to take, how to use?

"Almag-01" is a small device consisting of a chain of elastic coils. Feedback from patients who have used the device indicates that it is easy to use and does not require special skills.



device almag 01 reviews